



## **SciDNA™ Size Standards**

Cat# :SciS-500

### **Description:**

Debna BioGene's DNA Size Standard series are used in assigning sizes to DNA fragments on fluorescence-detecting instruments. The common applications include genotyping and DNA fragment analysis. Each of these standards consist of 15 DNA fragments, ranging in 35, 50, 75, 100, 139, 150, 160, 200, 250, 300, 340, 350, 400, 450, 490 and 500 bp. Size fragments are evenly distributed and can be used for very accurate size calling.

### **Application:**

SciDNA™ Size Standards can be used at the same setting as of ABI's GeneScan 500 Liz Size Standard.

**Recommended Loading:** 0.5 µl per well.

**Storage temperature:** 4 °C for up to 1 month; -20 °C for 6 months.

For Research Use Only. Not for use in diagnostic procedures.